



RLT-S22-WT



HOME THERAPY RED LIGHT WAIST & BACK

User Manual

Relax Muscles + Enhance Circulation

Accelerate Healing + Relieve Pain

Enhance Recovery + Circulation

before using this product read this user manual carefully and keep for reference

PACKAGE CONTENTS



1. 1x Red & Infrared LED Light Therapy waist & back(2 x 108 pcs Red lights & 108 pcs IR lights.)

2. Remote (with 20-minutes auto shut-off timer),We provide two light modes for the device (steady lights / flickering lights), please switch the light mode by pressing the ON/OFF button.

3. 1x AC Adapter(AC 100-240V) & Power Cord

4. 1x DC Cord


Model:	RLT-S22-WT
Material:	75% neoprene and 25% velvet
Effective area:	12 X 6inch
Power supply:	AC100-240 V or DC 12V DGYAO Portable Battery
Frequency:	50/60 Hz
Wavelength:	660nm+880nm
Number:	2 x 108pcs 660nm red light & 108 pcs 880nm near infrared
Light Beads Craft:	SMD
Temperature :	Max 104 °F
N.W:	200g
G.W:	600g

INTRODUCTION


Thank you for purchasing your DGYAO red & infrared light therapy!

Each DGYAO Therapy Device pads has a built-in infrared lights (880nm wavelength)+red lights (660nm wavelength).circulation, accelerates nerve regeneration, reduces pain, relaxes muscles, and more.

We recommend using your device at least once a day for 20 minutes periodically to enjoy the full effect.



GETTING STARTED



Follow the below instructions to get started and enjoy the benefits of your new device:

1. Connect the DGYAO Therapy Device pads to either one of the two DC ports on the remote control *Note: Each remote has two DC ports. The 2nd port is a spare port.
2. Connect the AC Adapter to the remote.
3. Connect power cord to the AC Adapter.
4. Plug the power cord to the power socket (the plug adaptor needed for your country will be provided)
5. Instructions for use (two different light modes)
 - 5.1 Steady Lights Mode : Turn on the lights by pressing the ON/OFF button on the remote control, the lights will shut off automatically after 20 minutes' treatment. You may turn it on again by pressing the ON/OFF button. Press down on the ON/OFF button for about 3 seconds to turn off the lights.
 - 5.2 Flickering Lights Mode: Turn on the lights by pressing the ON/OFF button on the remote control, then press the ON/OFF button again , the device will switch to flickering lights mode(pulse rate :10 HZ). The lights will shut off automatically after 20 minutes' treatment. Please note that when the two light modes are used alternately, the total time is 20 minutes.
6. For maximum benefit, We recommend using the light therapy device pads on naked skin directlyusing the light therapy Devices pads on naked skin directly.

RED LIGHT THERAPY WARRANTY

Limited One-Year Warranty

DGYAO carries a limited one-year warranty for the original purchaser against any technical defects. Warranty coverage applies to the original owner and to the original product only and is not transferable.

DGYAO promises that the device has one year warranty for any defects in material and workmanship from the date of purchase when the device is used and maintained according to the requirements outlined in the User Manual.


If your device stops functioning or working as intended during the one-year warranty period, please call DGYAO Customer Care at Skype +86-139-2922-1675 and one of our customer service team members will work with you to solve the issue and give you instructions on how to get your DGYAO LIGHT THERAPY repaired or replaced, free of charge.

This warranty does not cover normal wear of the device and does not apply to any device that has been tampered with or damaged due to misuse, abuse, or negligent handling.

Thanks for your cooperation in advance.

SPECIFIC WARRANTY POLICY

- 1-Year defects liability period.
- 60 days return and replacement warranty and 18-month free replacement for charger, controller and connectors. (if they are broken). Return postage beared by customers.
- Return requests as " I don't like it", "I don't want it", "the color doesn't look good" are not accepted.
- The situation of less than three light broken will not covered in the warranty policy because it doesn't affect usage effect.
- Keep the original outer package for successful return.



CAUSED BY
MANUAL LABOR



INTENSE
EXERCISE



SEDENTARY
AT WORK



THE ELDERLY

FREQUENTLY ASKED QUESTIONS

Q1: What's the wavelenath of the liahts?

A1: 660nm in red lights and 880nm in infrared lights.

Q2: Does it has both red and infrared lights?

A2: With 108pcs near infrared lights + 2 x 108pcs red lights.

Q3: Can I use this device in Europe?

A3: Yes, you only need a EU adaptor and it can be used in Europe.

Q4: Does it need to be plugged into the socket?

A4: Yes. There is also another model that is directly connected to the DGYAO power bank and can be used similarly(Sold separately).

Q5: Can multiple people share the device?

A5: Yes.The device can be used by different people and your pets too if you like.

Q6: What voltage can this device be used?

A6: The device can be used within AC110-240V.

Q7: Is it more effective on bare skin than over cloth?

A7: Yes, it is more effective on bare skin than over cloth.

Q8: Are there side effects?

A8: Don't worry about it, but please notice you may will have skin irritation if you 're allergic to the fibers.

Q9: How long will it take to experience relief of device?

A9: We recommend wearing it for 2-3 therapy sessions perday for the first month.20 minutes each time.Take arest and drink some water to relax yourself before next following use.Please be patient because it may take several weeks to experience optimal relief.

A10: Can people with heart disease or heart pacemaker use it?

A10: We recommend that you consult your primary care provider or physician when introducing Light Therapy.

- Contact your Doctor if you have any medical questions regarding use of this device.Contact your Doctor prior to using-taking any light.sensitive medications or medical questions.

- Light Sensitivity Test- the User Guide also recommends testing yourself for light sensitivity prioerto initial treatment.

Q11: What situations should I not use the therapy?

A11: Pain which may benefit from heat can be relieved by red and infrared light therapy.But under some situation,do'NOT use the light therapy:

1. High Fever.Infrared light may increase the temperature of your tissue and can worsen fever.

2. Tumour &tuberculosis.Improper use of infraredradia tion may lead to more severe tumors.
3. Patients fitted with heart pacemaker can NOT use light therapy.
4. Acute soft Tissue Injury. Bruises, sprains,bruises, falls or sport injury cause subcutaneous soft tissue torn bleeding or exudation.

For more information, please go to the blog of our official website: www.dgyao.net

Contact Us

- Consumers and clinicians call: +86-139-2922-1675 (Skype)- for all other questions or customer care needs, please email: dgyao@dgyao.net
- Customer Care Hours: 6am to 11am, 6pm-12pm Eastern Standard Time Monday – Sunday Except major holidays.
- Manufacturers' warranties may not apply in all cases, Please review the warranty carefully, and contact the manufacturer if you have any questions.

DONG GUAN TOP YAO INDUSTRY CO.,LTD

▶ Email: dgyao@dgyao.net

▶ Facebook: @dgyao

▶ Skype: +86 139 2922 1675

▶ Amazon Store: <https://www.amazon.com/dgyao>

▶ Official Website: <https://www.dgyao.net>



www.dgyao.net



FACEBOOK



TWITTER



YOUTUBE



LINKEDIN